

# RecoverAir 200

*Alternating Pressure  
Mattress Overlay*



## **Versatile, easy to use, cost effective pressure redistribution therapy**

For patients at high risk for skin breakdown, the RecoverAir 200 therapeutic mattress overlay provides alternating pressure every five minutes for lower interface pressure. All air cushions are constructed of the highest quality and most pliable coated nylon for low shear and increased patient comfort. A vapor permeable topcover supplies a low friction, low shear barrier and is easily washable. The RecoverAir 200 overlay is compatible with most facility mattresses.

**RecoverCare<sup>®</sup>**  
*Equipment. Solutions. Results!*

# RecoverAir 200 *Alternating Pressure Mattress Overlay*

## Features

- Alternating pressure every five minutes
- Inflation pressure calibrated to patient's weight
- Automatically adjusts to maintain optimal pressure redistribution
- Effective pressure redistribution throughout entire alternating cycle
- Hose quick release for rapid CPR deflation
- Compact, lightweight and made in the USA
- Available with optional RC Safety Perimeter system for reduced risk of falls
- Latex free



## Clinical Considerations

- Multiple stage 2, 3 or 4 pressure ulcers
- Contra-indicated for patients with unstable spine conditions

## Specifications

- Dimensions
  - Pump: 9" x 10" x 5" (5 lbs)
  - Overlay: (3.5 - 6") x 34" x 80"
- AP Cycle Time: 5 minutes
- Electrical rating: 120V/60Hz, 1A
- Agency Approvals: UL 544, C22.2 No. 125

### For Ordering Information

PHONE 888-750-7828 • FAX 866-750-7828

[www.recovercare.com](http://www.recovercare.com)

1920 Stanley Gault Parkway, Suite 100, Louisville, KY 40223

**RecoverCare**<sup>®</sup>  
Equipment. Solutions. Results!