

RecoverAir 2000

*Low Air Loss and
Alternating Pressure Overlay*



Progressive pressure redistribution therapy in a versatile overlay system

The RecoverAir 2000 uses both alternating pressure and low air loss to treat all stages of pressure ulcers. Progressive alternating pressure produces lower interface pressures to stimulate healing. Low Air Loss provides moisture and temperature control and improved patient comfort. The top cover is highly vapor permeable and oversized, thus reducing friction, shear and complications of incontinence. As a mattress overlay, the RecoverAir 2000 can be used with existing patient bed frames and mattresses.

RecoverCare[®]
Equipment. Solutions. Results!

RecoverAir 2000 *Low Air Loss and Alternating Pressure Overlay*

Features

- Provides effective low air loss to assist in moisture management
- Selectable alternating pressure with five minute cycle time
- Maximum Inflate for easy patient positioning
- Quilted top sheet is highly vapor permeable and fluid proof
- Overlay can be used with most facility bed frames and mattresses
- CPR quick release for rapid deflation
- Latex free

Clinical Considerations

- Multiple stage 2, 3 or 4 pressure ulcers
- Contra-indicated for patients with unstable spine conditions



Specifications

- Dimensions
 - Pump: 5" x 10" x 12" (9.5 lbs)
 - Mattress: 80" x 34" x 3.5" - 6.5"
- Weight limit: 350 lbs
- Electrical rating: 120V~, 60Hz, 5A/230V~, 50Hz, 2.5A
- Agency Approvals: UL 2601.1, C22.2; TUV, EN 60601-1

For Ordering Information

PHONE 888-750-7828 • FAX 866-750-7828

www.recovercare.com

1920 Stanley Gault Parkway, Suite 100, Louisville, KY 40223

RecoverCare[®]
Equipment. Solutions. Results!