

# RecoverTurn 3000

*Therapeutic Lateral  
Rotation Mattress  
Replacement System*



## **Progressive pulmonary and wound care therapy for critically ill patients**

The RecoverTurn 3000 combines a sophisticated lateral rotation system with advanced wound care therapy to provide treatment for the most critically ill patients. The easy to use, reliable, microprocessor driven control unit allows for customized rotation angles and times to maximize pulmonary outcomes. Therapeutic low air loss helps manage moisture and provides alternating pressure therapy to prevent and treat pressure ulcers. As a mattress replacement system, the RecoverTurn 3000 can fit most patient bed frames thus reducing costs associated with an entire full frame bed system.

**RecoverCare<sup>®</sup>**  
*Equipment. Solutions. Results!*

# RecoverTurn 3000 *Therapeutic Lateral Rotation Mattress Replacement System*

## Features

- Microprocessor designed for simplicity, reliability and ease of use
- Rotation angles are customized and set independently
- Rotation hold/dwell times are adjustable from 0 to 95 minutes
- Equipped with *Max Inflate* for easier patient egress or transfer
- Keypad may be locked to prevent interruption of therapy
- CPR quick release provides rapid deflation
- Provides uncompromised pressure redistribution throughout turning sequence
- *Pulsation* mode provides alternating pressure therapy
- Unique *Cradle* feature supports patient position
- Provides low air loss for moisture management
- Outer air manifolds help stabilize larger patients during the rotation cycle while providing additional protection from contact with side rails
- All air cushions are constructed of high quality and very pliable coated nylon for low shear and even lower interface pressures
- The quilted top sheet is highly vapor permeable and oversized providing a low friction, low shear barrier against complications due to moisture
- Convoluted foam base provides additional safety to your patient
- Latex free



## Clinical Considerations

- Acute intervention for pulmonary compromised patients
- All stages of pressure ulcers
- Contra-indicated for agitated patients or those with unstable spine conditions

## Specifications

- Dimensions
  - Pump: 6.5" x 10" x 12" (15.6 lbs)
  - Mattress: 80" x 34" x 10" (Larger mattress sizes available)
- Weight limit: 350 lbs (Bariatric sizes support up to 1000 lbs)
- Electrical rating: 115V/230V~,50Hz/60Hz, 5A/3A
- Agency Approvals: UL 2601.1, C22.2 No. 601.1

### For Ordering Information

PHONE 888-750-7828 • FAX 866-750-7828

[www.recovercare.com](http://www.recovercare.com)

1920 Stanley Gault Parkway, Suite 100, Louisville, KY 40223

**RecoverCare**<sup>®</sup>  
Equipment. Solutions. Results!