

RecoverAir 250

*Alternating Pressure
Mattress Replacement
System*



Active alternating pressure for cost-effective wound care treatment

The RecoverAir 250 therapeutic support surface provides alternating pressure every five minutes for lower interface pressure and promoted circulation. All air cushions are constructed of the highest quality and most pliable coated nylon for low shear and friction protection. The convoluted foam base provides additional safety to your patient, and the top sheet is highly vapor permeable and oversized to provide a low friction, low shear barrier against complications due to incontinence. Its modular design facilitates cleaning and interchangeability. As a mattress replacement system, the RecoverAir 250 can be used with existing facility bed frames.

RecoverCare[®]
Equipment. Solutions. Results!

RecoverAir 250 *Alternating Pressure Mattress Replacement System*

Features

- Active alternating pressure
- Fully adjustable to the patient's weight, extremely quiet
- Simple dial weight adjustment
- Programmed never to exceed 32mmHg inflation pressure
- CPR quick release for rapid deflation
- Compact, lightweight and made in the USA
- Latex free

Clinical Considerations

- Multiple stage 2, 3 or 4 pressure ulcers
- Contra-indicated for patients with unstable spine conditions



Specifications

- Dimensions
 - Pump: 9" x 10" x 5" (5 lbs)
 - Mattress: 10" x 34" x 80" (15 lbs)
- AP Cycle Time: 5 minutes
- Electrical rating: 120V/60Hz, 1A
- Agency Approvals: UL 544, C22.2 No. 125

For Ordering Information

PHONE 888-750-7828 • FAX 866-750-7828

www.recovercare.com

1920 Stanley Gault Parkway, Suite 100, Louisville, KY 40223

RecoverCare[®]
Equipment. Solutions. Results!