

RecoverAir 5000

*Zoned Low Air Loss
Mattress Replacement
System*



Multi-zone therapy for highly advanced wound treatment

The RecoverAir 5000 mattress replacement system offers advanced pressure redistribution and treatment for the most severe pressure ulcers. A sophisticated microprocessor provides multi-zoned pressure redistribution in which each zone can be adjusted independently for customized therapy. An optional low air loss button manages moisture. The lightweight RecoverAir 5000 mattress and control unit can be used with most patient bed frames to maximize your clinical and economic outcomes.

RecoverCare[®]
Equipment. Solutions. Results!

RecoverAir 5000 *Zoned Low Air Loss Mattress Replacement System*

Features

- Each pressure zone can be adjusted independently for customized therapy
- *Auto Profile* sets mattress pressure based on patient's height and weight
- Automatic transport mode maintains inflation during power outages or transport
- Microprocessor controlled for simplicity, reliability, and ease of use
- Keypad may be locked to prevent interruption of therapy
- Control unit is compact and lightweight
- Zone deflation is possible with a single keystroke
- Max inflate provides a firm surface for easy positioning
- CPR quick release for rapid deflation
- Each air cell provides more than eight inches of air support
- Air cushions are constructed of high quality and very pliable nylon for low shear and friction protection
- Convuluted foam base provides for additional safety to your patient
- Optional low air loss can be activated to manage moisture
- Quilted top sheet is highly vapor permeable and oversized providing a low friction, low shear barrier against complications due to moisture
- Latex free



Clinical Considerations

- Multiple stage 2, 3 or 4 pressure ulcers
- Contra-indicated for patients with unstable spine conditions

Specifications

- Dimensions
 - Pump: 6" x 10" x 12" (15.2 lbs)
 - Mattress: 80" x 34" x 10"
- Weight limit: 350 lbs
- Electrical rating: 115V/230V~,50Hz/60Hz, 5A/3A
- Agency Approvals: UL 2601.1, C22.2 No. 601.1

For Ordering Information

PHONE 888-750-7828 • FAX 866-750-7828

www.recovercare.com

1920 Stanley Gault Parkway, Suite 100, Louisville, KY 40223

RecoverCare[®]
Equipment. Solutions. Results!